

<b>TestCase Validation Report</b>	<b>04 01 2016, 13:03:09.198-04:00</b>
-----------------------------------	---------------------------------------

<b>Validation Type</b>	Context-Based
------------------------	---------------

<b>Result</b>	Failed - Not Supported
---------------	------------------------

<b>Comments</b>
TC comments

<b>Test Case</b>	Test Plan	TestPlan-Unit-Testing
	Test Group	
	Test Case	TestCase-Unit-Testing

<b>1 - TestStep Validation Report</b>	<b>04 01 2016, 13:03:09.176-04:00</b>
---------------------------------------	---------------------------------------

<b>Name</b>	Transmit the immunization report to the Immunization Registry
-------------	---

<b>Validation Type</b>	
------------------------	--

<b>Type</b>	SUT_INITIATOR
-------------	---------------

<b>Result</b>	Failed
---------------	--------

<b>Comments</b>
TestStep1 comments

<b>2 - TestStep Validation Report</b>	<b>04 01 2016, 13:03:09.177-04:00</b>
---------------------------------------	---------------------------------------

<b>Name</b>	Receive Acknowledgement
-------------	-------------------------

<b>Validation Type</b>	
------------------------	--

<b>Type</b>	TA_RESPONDER
-------------	--------------

<b>Result</b>	Failed - Not Supported
---------------	------------------------

<b>Comments</b>

TestStep2 comments

<b>Profile</b>	Name	Return an Acknowledgement
	Organization	NIST
	Type	Constrainable
	Profile Version	1.0.0
	Profile Date	20160219
	Standard	2.5.1

### Message

MSH|^~\&|NIST APP|NIST Test Iz Reg|Test EHR Application|X68|200906040000-0500||ACK^V04^ACK  
|NIST-IZ-ACK-001.01|T|2.5.1|||NE|NE||||Z23^CDCPHINVS

MSA|AA|NIST-IZ-001.01

### Failures interpretation

Category	Classification
RE_USAGE	WARNING
CONSTRAINT_FAILURE	ERROR
CODED_ELEMENT	ERROR
CONTENT_SPEC_ERROR	AFFIRMATIVE
INVALID_CONTENT	ERROR
CONTENT_SUCCESS	AFFIRMATIVE
BINDING_LOCATION	AFFIRMATIVE
X_USAGE	ERROR
W_USAGE	ERROR
CONSTRAINT_SUCCESS	AFFIRMATIVE
CONSTRAINT_SPEC_ERROR	AFFIRMATIVE
PREDICATE_SPEC_ERROR	AFFIRMATIVE
VALUE_SET	ERROR
VALUE_SET	ALERT
VALUE_SET	WARNING
EXTRA	ERROR
LENGTH	ERROR
PREDICATE_SUCCESS	AFFIRMATIVE
PREDICATE	ERROR
R_USAGE	ERROR
UNEXPECTED	ERROR
FORMAT	ERROR
CONTENT_FAILURE	ERROR
UNESCAPED_SEPARATOR	ERROR
CONFORMANCE_STATEMENT	ERROR
CARDINALITY	ERROR

Summary	
<b>1 Errors</b>	
<b>9 Alerts</b>	
<b>3 Warnings</b>	

Validation Errors		Count : 1
CONSTRAINT_FAILURE		Count : 2
1	Type :	CONSTRAINT_FAILURE
	Description:	IZ-TS_Z - The value of MSH-7.1 (Time) SHALL be formatted with YYYYMMDDHHMMSS+-ZZZZ.
	Location:	Line: 1 Column: 61 Path: MSH[1]-7[1].1
2	Type :	CONSTRAINT_FAILURE
	Description:	IZ-TS_Z - The value of MSH-7.1 (Time) SHALL be formatted with YYYYMMDDHHMMSS+-ZZZZ.
	Location:	Line: 1 Column: 61 Path: MSH[1]-7[1].1

Validation Alerts		Count : 9
VALUE_SET		Count : 18
1	Type :	VALUE_SET
	Description:	The value set HL70363_IZ has been excluded from the validation
	Location:	Line: 1 Column: 10 Path: MSH[1]-3[1]
2	Type :	VALUE_SET
	Description:	The value set HL70363_IZ has been excluded from the validation
	Location:	Line: 1 Column: 10 Path: MSH[1]-3[1].1
3	Type :	VALUE_SET
	Description:	The value set HL70362_IZ has been excluded from the validation
	Location:	Line: 1 Column: 19 Path: MSH[1]-4[1]
4	Type :	VALUE_SET
	Description:	The value set HL70363_IZ has been excluded from the validation
	Location:	Line: 1 Column: 19

		<b>Path: MSH[1]-4[1].1</b>
5	Type :	VALUE_SET
	Description:	The value set HL70361_IZ has been excluded from the validation
	Location:	<b>Line: 1</b>
		<b>Column: 36</b> <b>Path: MSH[1]-5[1]</b>
6	Type :	VALUE_SET
	Description:	The value set HL70363_IZ has been excluded from the validation
	Location:	<b>Line: 1</b>
		<b>Column: 36</b> <b>Path: MSH[1]-5[1].1</b>
7	Type :	VALUE_SET
	Description:	The value set HL70362_IZ has been excluded from the validation
	Location:	<b>Line: 1</b>
		<b>Column: 57</b> <b>Path: MSH[1]-6[1]</b>
8	Type :	VALUE_SET
	Description:	The value set HL70363_IZ has been excluded from the validation
	Location:	<b>Line: 1</b>
		<b>Column: 57</b> <b>Path: MSH[1]-6[1].1</b>
9	Type :	VALUE_SET
	Description:	The value set HL70363_IZ has been excluded from the validation
	Location:	<b>Line: 1</b>
		<b>Column: 135</b> <b>Path: MSH[1]-21[1].2</b>
10	Type :	VALUE_SET
	Description:	The value set HL70363_IZ has been excluded from the validation
	Location:	<b>Line: 1</b>
		<b>Column: 10</b> <b>Path: MSH[1]-3[1]</b>
11	Type :	VALUE_SET
	Description:	The value set HL70363_IZ has been excluded from the validation
	Location:	<b>Line: 1</b>
		<b>Column: 10</b> <b>Path: MSH[1]-3[1].1</b>
12	Type :	VALUE_SET
	Description:	The value set HL70362_IZ has been excluded from the validation
	Location:	<b>Line: 1</b>
		<b>Column: 19</b> <b>Path: MSH[1]-4[1]</b>

13	Type :	VALUE_SET
	Description:	The value set HL70363_IZ has been excluded from the validation
	Location:	<b>Line: 1</b> <b>Column: 19</b> <b>Path: MSH[1]-4[1].1</b>
14	Type :	VALUE_SET
	Description:	The value set HL70361_IZ has been excluded from the validation
	Location:	<b>Line: 1</b> <b>Column: 36</b> <b>Path: MSH[1]-5[1]</b>
15	Type :	VALUE_SET
	Description:	The value set HL70363_IZ has been excluded from the validation
	Location:	<b>Line: 1</b> <b>Column: 36</b> <b>Path: MSH[1]-5[1].1</b>
16	Type :	VALUE_SET
	Description:	The value set HL70362_IZ has been excluded from the validation
	Location:	<b>Line: 1</b> <b>Column: 57</b> <b>Path: MSH[1]-6[1]</b>
17	Type :	VALUE_SET
	Description:	The value set HL70363_IZ has been excluded from the validation
	Location:	<b>Line: 1</b> <b>Column: 57</b> <b>Path: MSH[1]-6[1].1</b>
18	Type :	VALUE_SET
	Description:	The value set HL70363_IZ has been excluded from the validation
	Location:	<b>Line: 1</b> <b>Column: 135</b> <b>Path: MSH[1]-21[1].2</b>

Validation Warnings		Count : 3
<b>RE_USAGE</b>		<b>Count : 6</b>
1	Type :	RE_USAGE
	Description:	Segment ERR (Error) is missing. Depending on the use case and data availability it may be appropriate to value this element.(Usage is RE, Required, but may be Empty)
	Location:	<b>Line: 1</b> <b>Column: 1</b> <b>Path: ERR[1]</b>
2	Type :	RE_USAGE
	Description:	Field MSH-23 (Receiving Responsible Organization) is missing. Depending

	Location:	on the use case and data availability it may be appropriate to value this element.(Usage is RE, Required, but may be Empty) <b>Line: 1</b> <b>Column: 1</b> <b>Path: MSH[1]-23[1]</b>
3	Type : Description:  Location:	RE_USAGE Field MSH-22 (Sending Responsible Organization) is missing. Depending on the use case and data availability it may be appropriate to value this element.(Usage is RE, Required, but may be Empty) <b>Line: 1</b> <b>Column: 1</b> <b>Path: MSH[1]-22[1]</b>
4	Type : Description:  Location:	RE_USAGE Segment ERR (Error) is missing. Depending on the use case and data availability it may be appropriate to value this element.(Usage is RE, Required, but may be Empty) <b>Line: 1</b> <b>Column: 1</b> <b>Path: ERR[1]</b>
5	Type : Description:  Location:	RE_USAGE Field MSH-23 (Receiving Responsible Organization) is missing. Depending on the use case and data availability it may be appropriate to value this element.(Usage is RE, Required, but may be Empty) <b>Line: 1</b> <b>Column: 1</b> <b>Path: MSH[1]-23[1]</b>
6	Type : Description:  Location:	RE_USAGE Field MSH-22 (Sending Responsible Organization) is missing. Depending on the use case and data availability it may be appropriate to value this element.(Usage is RE, Required, but may be Empty) <b>Line: 1</b> <b>Column: 1</b> <b>Path: MSH[1]-22[1]</b>

<b>3 - TestStep Validation Report</b>	<b>04 01 2016, 13:03:09.188-04:00</b>
---------------------------------------	---------------------------------------

<b>Name</b>	Record an adverse reaction
-------------	----------------------------

<b>Validation Type</b>	
------------------------	--

<b>Type</b>	TA_RESPONDER
-------------	--------------

<b>Result</b>	Passed - Notable Exception
---------------	----------------------------

Comments
TestStep3 comments

<b>Profile</b>	Name	Return an Acknowledgement
	Organization	NIST
	Type	Constrainable
	Profile Version	1.0.0
	Profile Date	20160219
	Standard	2.5.1

Message
MSH ^~\& NIST APP NIST Test Iz Reg Test EHR Application X68 200906040000-0500  ACK^V04^ACK  NIST-IZ-ACK-001.01 T 2.5.1   NE NE    Z23^CDCPHINVS  MSA AA NIST-IZ-001.01

Failures interpretation	
Category	Classification
RE_USAGE	WARNING
CONSTRAINT_FAILURE	ERROR
CODED_ELEMENT	ERROR
CONTENT_SPEC_ERROR	AFFIRMATIVE
INVALID_CONTENT	ERROR
CONTENT_SUCCESS	AFFIRMATIVE
BINDING_LOCATION	AFFIRMATIVE
X_USAGE	ERROR
W_USAGE	ERROR
CONSTRAINT_SUCCESS	AFFIRMATIVE
CONSTRAINT_SPEC_ERROR	AFFIRMATIVE
PREDICATE_SPEC_ERROR	AFFIRMATIVE
VALUE_SET	ERROR
VALUE_SET	ALERT
VALUE_SET	WARNING
EXTRA	ERROR
LENGTH	ERROR
PREDICATE_SUCCESS	AFFIRMATIVE
PREDICATE	ERROR
R_USAGE	ERROR
UNEXPECTED	ERROR
FORMAT	ERROR
CONTENT_FAILURE	ERROR

UNESCAPED_SEPARATOR	ERROR
CONFORMANCE_STATEMENT	ERROR
CARDINALITY	ERROR

Summary
<b>1 Errors</b>
<b>9 Alerts</b>
<b>3 Warnings</b>

Validation Errors	Count : 1
-------------------	-----------

Validation Alerts	Count : 9
-------------------	-----------

Validation Warnings	Count : 3
---------------------	-----------